



READY TO USE

EXHAUSTED
(TO BE REGENERATED)

TECHNICAL INFORMATION

Dessicant Material:	Disidry's Orange Indicating Silica Gel
Packaging Material:	Aluminium (Container) / Paper / Stickers
Aluminium Foil Thickness:	0,5 mm
External Diameter:	81 mm
Dry Keeper Thickness:	17 mm
Total Net Weight:	30 g
Internal Container Capacity:	50 gr (Silica Gel)
Holes:	196 (98 on each side)
Life:	Depending on use. Generally 2-4 months (once stickers removed).
Storage:	Do not leave outside, directly exposed to the sun, in open environment.

HOW TO USE

When the beads are dry, they appear orange in color, and when they have become saturated with moisture to about 15-20% by weight the beads will change to a light blue color. This type of indication enables users to visually confirm when the desiccant is needs to be replaced or regenerated.

Silica Gel's life is variable and it depends on the amount of moisture that absorbs. To verify if it's still working, check the color: the orange Silica Gel turns to blue to indicate its complete exhaustion.

Silica Gel can be regenerated by heating in the oven for 1-2 hours at 100°C (212°F).

Fits all cases.

Plaber s.r.l. - HPRC cases

Via Div Tridentina, 10
I-36061 Bassano del Grappa
Vicenza - Italy

Tel +39 0424 809150
Fax +39 0424 808548

HPRC USA

240 Anderson Avenue
07074 Moonachie
New Jersey - Usa

Tel +1 201 617 2010
Fax +1 201 617 2001

www.hprc.it
E: info@hprc.it

WARNINGS



**DO NOT
SWALLOW**



**DO NOT
WARM IN A
MICROWAVE**



**DO NOT TOUCH
WHEN IT'S HOT**

(AFTER REGENERATION
IN THE OVEN)



**DO NOT
OPEN**